



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hahnemühle pleated filter grade 597, 185 mm, pack of 100 pieces

**€55.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Due to the high purity recommended for the determination of substances and sample preparation for sensitive detection methods.

Very low mineral content (ash 0.08%).

Bosonders suitable for: Food control, beverage analysis, environmental monitoring.

## Additional Information

---

No.	HN-0134
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4649
EAN	4058072261238
Transport temperature	Room temperature
diameter	185 mm
Basis weight	85 g/m <sup>2</sup>
Heart mountain s/100 ml	155
Surface/ Structure	smooth
Variety Hahnemühle	597
Variety Munktell	292 41b
Variety Whatman	597
Type filter	Pleated filter
for type analysis	qualitative

