



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hahnemühle pleated filter grade 597, 240 mm, pack of 100 pieces

€68.40
plus VAT &
Shipping

Product Images



Description

Due to the high purity recommended for the determination of substances and sample preparation for sensitive detection methods.

Very low mineral content (ash 0.08%).

Bosonders suitable for: Food control, beverage analysis, environmental monitoring.

Additional Information

No.	HN-0135
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4650
EAN	4058072261245
Transport temperature	Room temperature
diameter	240 mm
Basis weight	85 g/m ²
Heart mountain s/100 ml	155
Surface/ Structure	smooth
Variety Hahnemühle	597
Variety Munktell	292 41b
Variety Whatman	597
Type filter	Pleated filter
for type analysis	qualitative

