



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Magnetic Filter Holder, 300 ml, for Laboratory Bottles

€151.40
plus VAT &
Shipping

Product Images



Description

One filter holder for different bottles.

The magnetic filter holder is suitable for membrane filters with a diameter of 47 mm. The user-friendly design allows one-hand operation and avoids possible damage to the membrane.

You can use the magnetic filter holder with laboratory bottles (GL 33/45). The connection to a vacuum is possible at any time due to the included olive. Due to the material PES the filter holder is autoclavable and has a high chemical resistance.

Continuous filtration is possible due to the filter lid.

Volume: 300 ml

Applications:

- Microbiological filtration tests
- General filtration of aqueous samples

Additional Information

No.	8-3001
Manufacturer (Brand)	neoLabLine
EAN	4058072086992
Transport temperature	Room temperature
Volume	300 ml
sterile	No
Material housing	Polyethersulfone (PES)
for DM hose inside 1	8 to 12 mm
Spare part	8-3004 8-3005 8-3006
Type connection	Hose olive
Type filter	Magnetic filter holder

