



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Magnetic filter holder, 500 ml, for aspirating bottles

€181.90
plus VAT &
Shipping

Product Images



Description

One filter holder for different bottles.

The magnetic filter holder is suitable for membrane filters with a diameter of 47 mm. The user-friendly design allows one-hand operation and avoids possible damage to the membrane.

You can use the magnetic filter holder with suction bottles. The connection to a vacuum is possible at any time due to the included olive. Due to the material PES the filter holder is autoclavable and has a high chemical resistance.

Continuous filtration is possible due to the filter lid.

Volume: 500 ml

Applications:

- Microbiological filtration tests
- General filtration of aqueous samples

Additional Information

No.	8-3002
Manufacturer (Brand)	neoLabLine
EAN	4058072087005
Transport temperature	Room temperature
Volume	500 ml
sterile	No
Spare part	8-3004 8-3005 8-3006 8-3007 8-3008
Material housing	Polyethersulfone (PES)
Type connection	Hose olive
Type filter	Magnetic filter holder
for DM hose inside 1	8 to 12 mm

