



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Magnetic Filter Holder, 500 ml, for Laboratory Bottles

**€181.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

One filter holder for different bottles.

The magnetic filter holder is suitable for membrane filters with a diameter of 47 mm. The user-friendly design allows one-hand operation and avoids possible damage to the membrane.

You can use the magnetic filter holder with laboratory bottles (GL 33/45). The connection to a vacuum is possible at any time due to the included olive. Due to the material PES the filter holder is autoclavable and has a high chemical resistance.

Continuous filtration is possible due to the filter lid.

**Volume:** 500 ml

Applications:

- Microbiological filtration tests
- General filtration of aqueous samples

## Additional Information

---

No.	8-3003
Manufacturer (Brand)	neoLabLine
EAN	4058072087012
Transport temperature	Room temperature
Volume	500 ml
sterile	No
Spare part	8-3004 8-3005 8-3006
Material housing	Polyethersulfone (PES)
Type connection	Hose olive
Type filter	Magnetic filter holder
for DM hose inside 1	8 to 12 mm

