



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Millipore vacuum filter unit w. funnel & bottle 150 ml, 12 pcs/pack

€1,512.00
plus VAT &
Shipping

Product Images



Description

Millipore Stericup™ filter units for ultra-fast and economical clarification or sterile filtration of tissue culture media, biological fluids or other aqueous solutions.

Millipore Express® PLUS membrane (PES) with support grid, firmly integrated into the filter unit. This allows for maximum flow rate, reduced foaming and protein denaturation.

Hopper capacity: 150, 250, 500 ml.

Grippable and stable container.

Suitable for work in sterile workbenches.

Additional Information

No.	MP-1027
Manufacturer (Brand)	Millipore
old neoLab article no.	7-9060
EAN	4058072335502
Transport temperature	Room temperature
Color	blue
Volume	150 ml
sterile	Yes
Material housing	Polystyrene (PS)
Material membrane	Polyethersulfone (PES)
Pore size	0.22 µm
for DM hose inside 1	8 - 12 mm
Type connection	Hose olive
Type filter	Vacuum filter units

