



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**qpore® filtration unit PES, 250/500 ml, 0.45 µm,  
diameter: 50 mm, 12 pcs/pack**

**€159.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

These sterile **qpore® filtration units** with polyethersulfone (PES) membrane have the fastest flow rates and very low protein binding, making them ideal for filtration of cell culture media. At the end of filtration, the sterile solution can be stored in the filtrate bottle until needed.

The scope of delivery of the complete unit includes the filter funnel with 0.45 µm PES membrane, a vacuum connection, a lid, the filtrate bottle and a cap.

Features:

- Green color coding
- Quick connector for tubing with inner Ø 8-12 mm
- Graduated filter top
- Sterile, pyrogen free, detergent free
- The filtration units are individually sterile packed to 12 pieces
- Packaging with cut open aid; thus can be opened without further aids

Fields of application:

- Sterile filtration of cell culture media
- Microbiological media
- Buffer solutions and other biological solutions

## Additional Information

---

No.	6-0055
Manufacturer (Brand)	qpore
EAN	4058072169596
Transport temperature	Room temperature
Color	green
Volume	500 ml
sterile	Yes
Material housing	Polypropylene (PP)
Material membrane	Polyethersulfone (PES)
Pore size	0.45 µm
Type connection	Hose olive
Type filter	Vacuum filter units
for DM hose inside 1	10

