



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## neolab® membrane filter PTFE, non-sterile, 0.45 µm, 25 mm Ø, 100 pcs./pack

Special Price

**€150.00** was

**€225.00**

**plus VAT &  
Shipping**

### Product Images

---

## Description

---

Membrane filters made of PTFE are permanently hydrophobic, ideal for filtration of air and gases. Excellent chemical compatibility, suitable for filtration of acids and aggressive solvents, for aeration and deaeration applications, phase separation, aerosol collection and gas filtration.

Steam sterilizable. Non-sterile.

Thickness: 65-100 µm depending on pore diameter

Flow rate: 0.45 µm = 0.3 l/min

## Additional Information

---

No.	8-1145
Manufacturer (Brand)	neoLab
EAN	4058072050351
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
sterile	No
DM outside	25 mm
TBST MAX	121 °C
Filter properties	with net print
Area diaphragm	4.9 cm <sup>2</sup>
Fluid behavior	hydrophobic
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.45 µm
Type filter	Membrane filter

