



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Millipore membrane filter 25 mm, 0.22 μ m,
white without net imprint, 100 pcs/pack**

**€145.00
plus VAT &
Shipping**

Product Images



Description

The standard for laboratory applications, e.g. for analytical applications, high purification and sterile filtration. The membranes are resistant to dilute acids and alkalis, autoclavable at 121°C.

Main applications :

Clarification filtration of sera, cell harvesting and washing, colony hybridization, contamination analysis and microbiological analysis.

PU: 100pcs./pack

Additional Information

No.	MP-0379
Manufacturer (Brand)	Millipore
old neoLab article no.	7-9030
EAN	4058072335373
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
sterile	No
DM outside	25 mm
TBST MAX	121 °C
Fluid behavior	hydrophilic
Material membrane	Cellulose mixed ester
Pore size	0.22 µm
Area diaphragm	4.9 cm ²
Type filter	Membrane filter

