



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Millipore membrane filter 47 mm, 0.45 µm,
black with mesh print, 100 pcs/pack**

**€294.00
plus VAT &
Shipping**

Product Images

Description

The standard for laboratory applications, e.g. for analytical applications, high purification and sterile filtration. The membranes are resistant to dilute acids and alkalis, autoclavable at 121°C.

Main applications :

Clarification filtration of sera, cell harvesting and washing, colony hybridization, contamination analysis and microbiological analysis.

PU: 100pcs./pack

Additional Information

No.	MP-0416
Manufacturer (Brand)	Millipore
old neoLab article no.	7-9040
EAN	4058072335397
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
sterile	No
DM outside	47 mm
TBST MAX	121 °C
Filter properties	with net print
Fluid behavior	hydrophilic
Material membrane	Cellulose mixed ester
Pore size	0.45 µm
Area diaphragm	17.3 cm ²
Type filter	Membrane filter

