



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Millipore membrane filter 25 mm, 0.45 µm,  
white with net imprint, 100 pcs/pack**

**€157.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

The standard for laboratory applications, e.g. for analytical applications, high purification and sterile filtration. The membranes are resistant to dilute acids and alkalis, autoclavable at 121°C.

**Main applications :**

Clarification filtration of sera, cell harvesting and washing, colony hybridization, contamination analysis and microbiological analysis.

**PU:** 100pcs./pack

## Additional Information

---

No.	MP-0431
Manufacturer (Brand)	Millipore
old neoLab article no.	7-9034
EAN	4058072335403
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
sterile	No
DM outside	25 mm
TBST MAX	121 °C
Filter properties	with net print
Area diaphragm	4.9 cm <sup>2</sup>
Fluid behavior	hydrophilic
Material membrane	Cellulose mixed ester
Pore size	0.45 µm
Type filter	Membrane filter

