



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Ahlstrom ReliaDisc Membrane Filter, Cellulose Nitrate (CN), white, without mesh, non-sterile

**€233.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Cellulose nitrate membranes are the standard in analytical and microbiological filtration. The membranes are well suited for microbiological analysis of water, food and beverages in the pH range 4–8. Maximum thermal stability: 130°C.

## Additional Information

---

No.	8-7008
Manufacturer (Brand)	Ahlstrom
EAN	4058072050986
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
sterile	No
DM outside	47 mm
TBST MAX	130 °C
Fluid behavior	hydrophilic
Material membrane	Cellulose nitrate
Pore size	1.20 µm
Flow rate water/isopropanol mixture per cm <sup>2</sup>	320 ml/min
Type filter	Membrane filter

