



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Ahlstrom ReliaDisc Membrane Filter, Cellulose Nitrate (CN), white, without mesh, non-sterile

€150.00
plus VAT &
Shipping

Product Images



Description

Cellulose nitrate membranes are the standard in analytical and microbiological filtration. The membranes are well suited for microbiological analysis of water, food and beverages in the pH range 4–8. Maximum thermal stability: 130°C.

Additional Information

No.	8-7011
Manufacturer (Brand)	Ahlstrom
EAN	4058072051013
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
sterile	No
DM outside	50 mm
TBST MAX	130 °C
Fluid behavior	hydrophilic
Material membrane	Cellulose nitrate
Pore size	0.45 µm
Flow rate water/isopropanol mixture per cm ²	70 ml/min
Type filter	Membrane filter

