



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Ahlstrom ReliaDisc membrane filter, Teflon (PTFE), non-sterile

**€570.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

ReliaDisc membrane filters made of Teflon (PTFE) are very well suited for particle filtration of air and gas. The permanently hydrophobic filter material has extreme resistance to aggressive solvents and acids (pH 1-14). Maximum thermal stability: 200°C. The hydrophobic filter must be moistened with ethanol before filtration of aqueous samples.

## Additional Information

---

No.	8-7015
Manufacturer (Brand)	Ahlstrom
EAN	4058072051051
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
sterile	No
DM outside	47 mm
TBST MAX	200 °C
Filter properties	with special chemical resistance
Fluid behavior	hydrophobic
Material membrane	Polytetrafluoroethylene (PTFE. Teflon®)
Pore size	0.45 µm
Flow rate isopropanol per cm <sup>2</sup>	0.2 ml/min
Type filter	Membrane filter

