



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Ahlstrom ReliaDisc membrane filter, polyamide (PA), non-sterile

€201.40
plus VAT &
Shipping

Product Images



Description

ReliaDisc membrane filters made of polyamide (PA)/nylon are widely used for filtration of aqueous (pH 3-14) as well as biological solutions. They are well suited for clear filtration of buffers and nutrient media with a low number of extractable components. The membrane filters are hydrophilic and have a maximum thermal stability of 100°C.

Additional Information

No.	8-7013
Manufacturer (Brand)	Ahlstrom
EAN	4058072051037
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
sterile	No
DM outside	47 mm
TBST MAX	100 °C
Flow rate water/isopropanol mixture per cm ²	35 ml/min
Fluid behavior	hydrophilic
Material membrane	Polyamide 6.6 (PA 6.6. Nylon®)
Pore size	0.45 µm
Type filter	Membrane filter

