



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Ahlstrom ReliaDisc membrane filter, regenerated cellulose (RC), non-sterile

€298.50
plus VAT &
Shipping

Product Images



Description

ReliaDisc regenerated cellulose (RC) membrane filters are hydrophilic and resistant to many solvents, excellent for removing particles from solvents, pH range 3-12. Typical application: sample preparation for HPLC. Maximum thermal stability: 200°C.

Additional Information

No.	8-7017
Manufacturer (Brand)	Ahlstrom
EAN	4058072051075
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
sterile	No
DM outside	47 mm
TBST MAX	200 °C
Fluid behavior	hydrophilic
Material membrane	Cellulose. regenerated (RC)
Pore size	0.45 µm
Flow rate water/isopropanol mixture per cm ²	30 ml/min
Type filter	Membrane filter

