



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**qpore® membrane filter, PTFE, non-sterile,  
hydrophilic, 0.22 µm, Ø 47 mm**

**€99.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

The hydrophilic qpore membrane filters made of polytetrafluoroethylene PTFE are ideally suited for the filtration of dipolar media where high chemical resistance is required. The PTFE membranes are biologically and chemically inert and exhibit high flux rates.

## Additional Information

---

|                       |                                |
|-----------------------|--------------------------------|
| No.                   | 6-0086                         |
| Manufacturer (Brand)  | qpore                          |
| EAN                   | 4058072173050                  |
| Transport temperature | Room temperature               |
| diameter              | 47 mm                          |
| sterile               | No                             |
| Fluid behavior        | hydrophilic                    |
| Material membrane     | Polytetrafluoroethylene (PTFE) |
| Pore size             | 0.22 µm                        |

