



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**qpore® membrane filter, PTFE, non-sterile,
hydrophilic, 0.45 µm, Ø 47 mm**

**€99.00
plus VAT &
Shipping**

Product Images



Description

The hydrophilic qpore membrane filters made of polytetrafluoroethylene PTFE are ideally suited for the filtration of dipolar media where high chemical resistance is required. The PTFE membranes are biologically and chemically inert and exhibit high flux rates.

Additional Information

No.	6-0087
Manufacturer (Brand)	qpore
EAN	4058072171667
Transport temperature	Room temperature
diameter	47 mm
sterile	No
Fluid behavior	hydrophilic
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.45 µm

