



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**qpore® membrane filter, PVDF, non-sterile,
hydrophilic, 0.45 µm, Ø 47 mm**

**€99.00
plus VAT &
Shipping**

Product Images



Description

The hydrophilic qpore membrane filters made of polyvinylidene fluoride PVDF bind far less protein than membrane filters made of nylon, PTFE or cellulose mixed ester. They have a very high flux rate compared to other membrane materials. Due to the low protein and low drug binding properties, membrane filters made of PVDF are excellent for filtration in the life science sector. These membrane filters are also excellent for clear filtration of air and other gases.

Additional Information

No.	6-0091
Manufacturer (Brand)	qpore
EAN	4058072173470
Transport temperature	Room temperature
diameter	47 mm
sterile	No
Fluid behavior	hydrophilic
Material membrane	Polyvinylidene fluoride (PVDF)
Pore size	0.45 µm

