



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hahnemühle round filter grade 400, 90 mm Ø, 100 pcs/pack

**€7.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Round filters for general laboratory work are recommended for substance identification, clarification of liquids and sample preparation. The filters are made of highly refined cellulose. Smooth surface, medium fast filtering.

Applications:

Sample preparations

Clarification of: Alcohols, essences, vinegar, essential oils, extracts/ gelatin, glycerin, paints, varnishes, hair tonics, perfumes, tinctures/ beer wort, spirits, syrup

Use as a protective sheet in filter presses

## Additional Information

---

No.	HN-0195
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4608
EAN	4058072261283
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	90 mm
Color	White
Material	Paper
sterile	No
Basis weight	61 g/m <sup>2</sup>
for type analysis	quantitative
Filtration speed	medium fast
Surface/ Structure	smooth
Variety Hahnemühle	400
Variety Munktell	4b
Variety Whatman	93
Type filter	Round filter

