



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hahнемühle round filter grade 400, 110 mm Ø, 100 pcs/pack

€8.80
plus VAT &
Shipping

Product Images



Description

Round filters for general laboratory work are recommended for substance identification, clarification of liquids and sample preparation. The filters are made of highly refined cellulose. Smooth surface, medium fast filtering.

Applications:

Sample preparations

Clarification of: Alcohols, essences, vinegar, essential oils, extracts/ gelatin, glycerin, paints, varnishes, hair tonics, perfumes, tinctures/ beer wort, spirits, syrup

Use as a protective sheet in filter presses

Additional Information

No.	HN-0196
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4609
EAN	4058072261290
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	110 mm
Color	White
Material	Paper
sterile	No
Basis weight	61 g/m ²
for type analysis	quantitative
Filtration speed	medium fast
Surface/ Structure	smooth
Variety Hahnemühle	400
Variety Munktell	4b
Variety Whatman	93
Type filter	Round filter

