




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Ahlstrom round filter qualitative, type 3 hw, 90
mm Ø, 100 pcs./pack**

**€8.84
plus VAT &
Shipping**

Product Images



Description

Munktell type 3hw,100pcs/pack

Additional Information

No.	E-1405
Manufacturer (Brand)	Ahlstrom
EAN	4058072064525
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	90 mm
Color	White
Material	Paper
sterile	No
Basis weight	65 g/m ²
for type analysis	qualitative
Filtration speed	20 sec/10 ml
Surface/ Structure	smooth
Variety Hahnemühle	2555
Variety Munktell	3 hw 15 1002
Variety Whatman	student grade
Type filter	Round filter

