



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hahnemühle round filter 589/6, 150 mm Ø, 100 pcs/pack

€85.30
plus VAT &
Shipping

Product Images



Description

Ash-free round filter papers for quantitative analysis.

The round filters are made of highly refined cotton linters.

Cellulose content over 95%, therefore high stability and durability.

Acid washed and neutralized with water. Free of minerals and metal ions, making them ideally suited for the determination of trace elements.

For analytical applications, routine quantitative and gravimetric procedures in food control, beverage analysis and environmental monitoring.

Applications:

For fine crystalline precipitates.;hot precipitates.

Additional Information

No.	HN-0428
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4641
EAN	4066292542495
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	150 mm
Color	green
Material	Paper
sterile	No
Basis weight	74 g/m ²
for type analysis	quantitative
Filtration speed	slow
Surface/ Structure	smooth
Variety Hahnemühle	589/6
Variety Munktell	390
Variety Whatman	44
Type filter	Round filter

