



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hahnemühle round filter 589/6, 185 mm Ø, 100 pcs/pack

**€85.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Ash-free round filter papers for quantitative analysis.

The round filters are made of highly refined cotton linters.

Cellulose content over 95%, therefore high stability and durability.

Acid washed and neutralized with water. Free of minerals and metal ions, making them ideally suited for the determination of trace elements.

For analytical applications, routine quantitative and gravimetric procedures in food control, beverage analysis and environmental monitoring.

Applications:

For fine crystalline precipitates.;hot precipitates.

## Additional Information

---

No.	HN-0429
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4642
EAN	4066292542501
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	185 mm
Color	green
Material	Paper
sterile	No
Basis weight	74 g/m <sup>2</sup>
for type analysis	quantitative
Filtration speed	slow
Surface/ Structure	smooth
Variety Hahnemühle	589/6
Variety Munktell	390
Variety Whatman	44
Type filter	Round filter

