



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Round Filter quantitative, Type 389, 110 mm Ø, 100 pcs./pack

€31.70
plus VAT &
Shipping

Product Images



Description

Munktell, ash-free
Paper weight: 80 g/sqm,
Filtering time: FS * 25
White point

Additional Information

No.	E-1745
Manufacturer (Brand)	neolab
EAN	4058072066888
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	110 mm
Color	White
Material	Paper
sterile	No
Basis weight	84 g/m ²
for type analysis	quantitative
Filtration speed	20 sec/10 ml
Surface/ Structure	smooth
Variety Hahnemühle	589/2
Variety Munktell	389
Variety Whatman	40 43
Type filter	Round filter

