



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Socorex protection filter made of PP fibers for  
macropipette 10 mL, Calibra® digital 832,  
Acura® manual 835 and electro 935/936**

**€47.32  
plus VAT &  
Shipping**

## Product Images

## Description

---

# Protective filters for Socorex pipettes

For protection against liquid ingress or contamination.

The filters are replaceable.

For additional safety against any liquid ingress and contamination risk.

Non-sterile, non-autoclavable.

## Additional Information

---

No.	SR-0413
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7870
EAN	4058072369262
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
suitable for	Single-channel pipettes Multi-channel pipettes
Fluid behavior	hydrophilic
Material membrane	PP fibers
Type filter	Safety filter

