



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Sartorius safety filter PP/PTFE 0.45 µm, 50 mm Ø, 99.97% recovery.

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Protection against harmful contamination on valves and vacuum pump parts!

High efficiency, effective retention of bacteria, low impact on performance. Protects air from contamination and oil aerosols.

Technical data

Pore size: 0.45 µm

Max. operating temperature: 134°C

Accessories for suction unit (Order No. 2-9335/36 and 2-9339)

Additional Information

No.	ST-3674
Manufacturer (Brand)	Sartorius
old neoLab article no.	3-1905
EAN	4058072371180
Transport temperature	Room temperature
DM outside	64 mm
TBST MAX	134 °C
MAX operating pressure	3 bar abs.
Material housing	Polypropylene (PP)
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.45 µm
Type Connection Output	Hose olive
Type Connection Input	Hose olive
for DM hose inside 1	6.4 to 12.7 mm
Area diaphragm	20 cm ²
Air flow rate	8.5 l/min (at 0.1 bar)
Type filter	Safety filter

