



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## qpore® Noliquids safety filter PTFE, 0.45 µm, Ø 50 mm

**€14.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**qpore®** Noliqids safety filter prevents contamination of room air with harmful aerosols and protects vacuum pumps from liquid penetration. Also suitable for sterile ventilation of bioreactors.

The filters can be autoclaved several times.

This keeps the costs low !

## Additional Information

---

No.	6-0036
Manufacturer (Brand)	qpore
EAN	4058072169640
Transport temperature	Room temperature
autoclavable at 121°C	Yes
DM outside	53 mm
TBST MAX	100 °C
MAX operating pressure	3.5 bar abs.
Material housing	Polypropylene (PP)
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.45 µm
Type Connection Output	Hose olive
Type Connection Input	Hose olive
for DM hose inside 1	8 - 12 mm
Area diaphragm	19.62 cm <sup>2</sup>
Air flow rate	12 l/min (at 0.1 bar)
Type filter	Safety filter

