





Macherey & Nagel CHROMAFIL PTFE-20/25

Product Images



Description

CHROMAFIL disposable filter PTFE-20/25 Membrane material: polytetrafluoroethylene (PTFE) Pore diameter: 0.2 µm, membrane diameter: 25 mm Color code top: yellow / bottom: colorless former product designation CHROMAFIL O-20/25 Pack of 100 pcs.

Additional Information

No.	ML-3076
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
EAN	4058072932992
Gross weight	0.362
Manufacturer EAN	4046681074462
Storage temperature	Raumtemperatur
Net weight	0.36200 kg
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.11300 m
sterile	No
Packing height	0.15000 m
Packing depth	0.11300 m
Packaging volume	0.00191500 ccm
Customs tariff number	84212980
Material membrane	Polytetrafluorethylen (PTFE)
Pore size	0.20 µm

