



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® filtration attachments sterile pore size  
0.2 µm, 28 mm Ø, CA membrane, 50 pcs./pack**

**€156.90  
plus VAT &  
Shipping**

## Product Images



## Description

---

Ready-to-use filter units.

Only virgin raw material is used for the ready-to-use filter units in order to exclude interaction with the sample to be filtered when used properly. Intended for laboratory use only.

**Diameter:** 28 mm

Individually sterile packed.

Connection: Inlet Luer lock inside, outlet Luer cone.

**PU:** 50 pcs./pack

## Additional Information

No.	E-1726
Manufacturer (Brand)	neoLab
EAN	4058072066772
Transport temperature	Room temperature
Color	blue
sterile	Yes
suitable for	Syringes
DM outside	28 mm
TBST MAX	45 °C
Filter properties	individually packed
MAX operating pressure	5 bar abs.
Fluid behavior	hydrophobic
Material housing	Polymethylmetacrylate (PMMA. PXM®) (acrylic glass)
Material membrane	Cellulose acetate (CA)
Pore size	0.20 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	5.3 cm <sup>2</sup>
Dead volume	0.1 ml
Type filter	Syringe pre-filter
for medium	Organic solutions

