



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® HPLC syringe filter, glass fiber prefilter
30 mm, 0.45 µm, PTFE, yellow, 100 pcs./pack**

**€244.40
plus VAT &
Shipping**

Product Images



Description

HPLC syringe filter with glass fiber pre-filter for pre-filtration of the coarsest particles to prevent membrane clogging. With additional edge overmolding for better pressure resistance of the filter unit and improved fixation of the membranes in the filter housing.

The imprint, together with the color coding, excludes any risk of confusion.

Effective filter area 4.91 cm²

Operating pressure 620 kPa (=6.2 bar/ 90 psi)

Inlet: Luer lock female

Outlet: Luer male

Characteristics

PTFE: hydrophobic, highly chemically resistant

Main applications

PTFE: organic solutions, aggressive acids and alkalis, cryogenic fluids

PU: 100pcs/pack

Additional Information

No.	7-8801
Manufacturer (Brand)	neoLab
EAN	4058072042677
Transport temperature	Room temperature
Color	yellow
sterile	No
suitable for	Syringes
DM outside	30 mm
TBST MAX	121 °C
Filter properties	with special chemical resistance with glass fiber prefilter suitable for HPLC
MAX operating pressure	6.2 bar abs.
Fluid behavior	hydrophobic
Material housing	Polypropylene (PP)
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.45 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	4.8 cm ²
Dead volume	0.14 ml
Type filter	Syringe pre-filter
for medium	Organic solutions

