



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® HPLC syringe filter, glass fiber prefilter 30 mm, 0.2 µm, RC, gray, 100 pcs./pack

€425.90
plus VAT &
Shipping

Product Images



Description

HPLC syringe filter with glass fiber pre-filter for pre-filtration of the coarsest particles to prevent membrane clogging. With additional edge overmolding for better pressure resistance of the filter unit and improved fixation of the membranes in the filter housing.

The imprint, together with the color coding, excludes any risk of confusion.

Effective filter area 4.91 cm²

Operating pressure 620 kPa (=6.2 bar/ 90 psi)

Inlet: Luer lock female

Outlet: Luer male

Characteristics

Regenerated cellulose (RC): hydrophilic, good wetting properties, very wet strength

Main applications

RC: aqueous and organic media (no acids/alkalis)

PU: 100pcs/pack

Additional Information

No.	7-8802
Manufacturer (Brand)	neoLab
EAN	4058072042684
Transport temperature	Room temperature
Color	gray
sterile	No
suitable for	Syringes
Filter properties	with special chemical resistance with glass fiber prefilter suitable for HPLC
MAX operating pressure	6.2 bar abs.
DM outside	30 mm
Area diaphragm	4.8 cm ²
Fluid behavior	hydrophilic
Material housing	Polypropylene (PP)
Material membrane	Cellulose. regenerated (RC)
Pore size	0.20 µm
TBST MAX	121 °C
Dead volume	0.14 ml
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Organic solutions Aqueous solutions

