



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® HPLC syringe filter, glass fiber prefilter  
30 mm, 0.45 µm, nylon, green, 100 pcs./pack**

**€247.10  
plus VAT &  
Shipping**

## Product Images



## Description

---

HPLC syringe filter with glass fiber pre-filter for pre-filtration of the coarsest particles to prevent membrane clogging. With additional edge overmolding for better pressure resistance of the filter unit and improved fixation of the membranes in the filter housing.

The imprint, together with the color coding, excludes any risk of confusion.

Effective filter area 4.91 cm<sup>2</sup>

Operating pressure 620 kPa (=6.2 bar/ 90 psi)

Inlet: Luer lock female

Outlet: Luer male

Characteristics

Nylon (PA): hydrophilic, very resistant to wet and dry conditions

Main applications

PA: aqueous and organic media (pharmaceutical applications)

**PU** : 100pcs/pack

## Additional Information

No.	7-8805
Manufacturer (Brand)	neoLab
EAN	4058072042714
Transport temperature	Room temperature
Color	green
sterile	No
suitable for	Syringes
DM outside	30 mm
TBST MAX	121 °C
Filter properties	with special chemical resistance with glass fiber prefilter suitable for HPLC
MAX operating pressure	6.2 bar abs.
Area diaphragm	4.8 cm <sup>2</sup>
Fluid behavior	hydrophilic
Material housing	Polypropylene (PP)
Material membrane	Polyamide 6.6 (PA 6.6. Nylon®)
Pore size	0.45 µm
Dead volume	0.14 ml
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Organic solutions Aqueous solutions

