



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Millipore syringe filter 0.22 μm , \varnothing 33 mm, 50 pcs/pack

€231.00
plus VAT &
Shipping

Product Images



Description

Order No. 7-9000

For sterile filtration of aqueous solutions.

Sample volume: < 100 ml

Packaging unit: 50 pcs./pack

Additional Information

No.	MP-1138
Manufacturer (Brand)	Millipore
old neoLab article no.	7-9000
EAN	4058072335564
Transport temperature	Room temperature
Color	blue
sterile	Yes
suitable for	Syringes
DM outside	33 mm
TBST MAX	45 °C
Filter properties	individually packed
MAX operating pressure	10 bar abs.
Area diaphragm	5.3 cm ²
Fluid behavior	hydrophilic
Material housing	Polymethylmetacrylate (PMMA, PXM®) (acrylic glass)
Material membrane	Cellulose mixed ester
Pore size	0.22 µm
Dead volume	0.1 ml
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Aqueous solutions

