



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Millipore syringe filter 0.45 $\mu\text{m}$ , $\varnothing$ 33 mm, 50 pcs/pack

**€236.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Order No. 7-9001**

For particle removal from aqueous solutions.

Sample volume: < 100 ml

**PU:** 50 pcs./pack

## Additional Information

No.	MP-1163
Manufacturer (Brand)	Millipore
old neoLab article no.	7-9001
EAN	4058072335595
Transport temperature	Room temperature
Color	blue
sterile	Yes
suitable for	Syringes
DM outside	33 mm
TBST MAX	45 °C
Filter properties	individually packed
MAX operating pressure	10 bar abs.
Area diaphragm	5.3 cm <sup>2</sup>
Fluid behavior	hydrophilic
Material housing	Polymethylmetacrylate (PMMA, PXM®) (acrylic glass)
Material membrane	Cellulose mixed ester
Pore size	0.45 µm
Dead volume	0.1 ml
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Aqueous solutions

