



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Millipore syringe filter 0.45 µm, Durapore, 33 mm Ø, 50 pcs/pack

€276.00
plus VAT &
Shipping

Product Images



Description

Order No. 7-9003

For particle removal from aqueous/organic solutions. Membrane with extremely low protein binding.
Sample volume: < 100 ml

PU: 50 pcs./pack

Additional Information

No.	MP-1193
Manufacturer (Brand)	Millipore
old neoLab article no.	7-9003
EAN	4058072335601
Transport temperature	Room temperature
Color	yellow
sterile	Yes
suitable for	Syringes
Filter properties	with special chemical resistance individually packed
MAX operating pressure	10 bar abs.
DM outside	33 mm
Area diaphragm	5.3 cm ²
Fluid behavior	hydrophilic
Material housing	Polymethylmetacrylate (PMMA. PXM®) (acrylic glass)
Material membrane	Polyvinylidene fluoride (PVDF)
Pore size	0.45 µm
TBST MAX	45 °C
Dead volume	0.1 ml
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Aqueous solutions Organic solutions
special properties	with extremely low protein binding

