



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Ahlstrom ReliaPrep syringe filter,  
polyethersulfone, D: 28 mm, sterile, pore size 0.2  
µm, Luer Lock, 50 pcs./pc**

**€196.10  
plus VAT &  
Shipping**

## Product Images



## Description

For sample preparation of aqueous and organic solutions, particle removal, high purification and pre-filtration, and sterile filtration. Top quality with high quality constructed housings. Polyethersulfone (PES) membrane material, sterile (ethylene oxide) and non-sterile versions. Polyethersulfone (PES): hydrophilic membrane with good chemical resistance and very high flow rates with very low protein binding. Housing material: methyl acrylate-butadiene-styrene polymer (MBS) Dead volume: 100 µl Filter area: 5.3 cm<sup>2</sup>

## Additional Information

No.	8-7063
Manufacturer (Brand)	Ahlstrom
EAN	4058072051532
Transport temperature	Room temperature
Color	transparent (clear)
sterile	Yes
suitable for	Syringes
DM outside	28 mm
Filter properties	with special chemical resistance individually packed
Area diaphragm	5.3 cm <sup>2</sup>
Fluid behavior	hydrophilic
Material housing	Acrylonitrile butadiene styrene (ABS)
Material membrane	Polyethersulfone (PES)
Pore size	0.22 µm
Dead volume	0.1 ml
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Liquids Organic solutions Aqueous solutions

