



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Ahlstrom ReliaPrep syringe filter, CA, D: 28 mm,  
sterile, pore size 0.2 µm, Luer Lock, 50 pcs./pc**

**€196.10  
plus VAT &  
Shipping**

## Product Images



## Description

For sample preparation of aqueous and organic solutions, particle removal, high purification and pre-filtration, and sterile filtration. Top quality with high quality constructed housings. Cellulose acetate membrane material, sterile (ethylene oxide) and non-sterile versions. Cellulose Acetate (CA): hydrophilic membrane with low protein binding, well suited for sterile and clear filtration of aqueous solutions, nutrient media and buffers. Housing material: methyl acrylate butadiene styrene polymer (MBS) Dead volume: 100 µl Filter area: 5.3 cm<sup>2</sup>

## Additional Information

No.	8-7043
Manufacturer (Brand)	Ahlstrom
EAN	4058072051334
Transport temperature	Room temperature
Color	transparent (clear)
sterile	Yes
suitable for	Syringes
DM outside	28 mm
Filter properties	with special chemical resistance individually packed
Fluid behavior	hydrophilic
Material housing	Acrylonitrile butadiene styrene (ABS)
Material membrane	Cellulose acetate (CA)
Pore size	0.22 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	5.3 cm <sup>2</sup>
Dead volume	0.1 ml
Type filter	Syringe pre-filter
for medium	Liquids Organic solutions Aqueous solutions

