



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Ahlstrom ReliaPrep syringe filter, glass fiber, D:
28 mm, non-sterile, pore size 0.7 µm, Luer Lock,
500 pcs./pc**

**€1,156.90
plus VAT &
Shipping**

Product Images



Description

For sample preparation of aqueous and organic solutions, particle removal, high purification and pre-filtration, and sterile filtration. Top quality with high quality constructed housings. Membrane material made of glass fiber. Glass fiber (GF): especially well suited for prefiltration of heavily loaded samples, high particle holding capacity. Housing material: methyl acrylate butadiene styrene polymer (MBS) Dead volume: 230 µl Filter area: 5.3 cm².

Additional Information

No.	8-7053
Manufacturer (Brand)	Ahlstrom
EAN	4058072051433
Transport temperature	Room temperature
Color	transparent (clear)
sterile	No
suitable for	Syringes
DM outside	28 mm
Filter properties	with special chemical resistance
Fluid behavior	hydrophilic
Material housing	Acrylonitrile butadiene styrene (ABS)
Material membrane	Fiberglass
Pore size	0.70 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	5.3 cm ²
Dead volume	0.23 ml
Type filter	Syringe pre-filter
for medium	Liquids Organic solutions Aqueous solutions

