



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Ahlstrom ReliaPrep syringe filter, nylon, D: 25 mm, sterile , pore size 0.2  $\mu$ m, Luer Slip, 50 pcs/pc**

**€181.80  
plus VAT &  
Shipping**

## Product Images



## Description

For sample preparation of aqueous and organic solutions, particle removal, high purification and pre-filtration, and sterile filtration. Top quality with high quality constructed housings. Nylon (NY) membrane material, sterile (ethylene oxide) and non-sterile versions. Nylon (NY): hydrophilic membrane with versatile application range with low content of extractable substances. Housing material: polypropylene (PP) Dead volume: 10 µl (15 mm), 100 µl (25 mm) Filter area: 1.7 cm<sup>2</sup> (15 mm) 4.8 cm<sup>2</sup> (25 mm).

## Additional Information

No.	8-7058
Manufacturer (Brand)	Ahlstrom
EAN	4058072051488
Transport temperature	Room temperature
Color	transparent (clear)
sterile	Yes
suitable for	Syringes
DM outside	25 mm
Filter properties	with special chemical resistance individually packed
Fluid behavior	hydrophilic
Material housing	Polypropylene (PP)
Material membrane	Polyamide 6.6 (PA 6.6. Nylon®)
Pore size	0.22 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer cone female
Area diaphragm	4.8 cm <sup>2</sup>
Dead volume	0.1 ml
Type filter	Syringe pre-filter
for medium	Liquids Organic solutions Aqueous solutions

