



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Ahlstrom ReliaPrep syringe filter, Teflon (PTFE)  
D: 4 mm, non-sterile, pore size 0.2  $\mu$ m, Luer Slip,  
500 pcs./pc**

**€711.80  
plus VAT &  
Shipping**

## Product Images



## Description

---

For sample preparation of aqueous and organic solutions, particle removal, high purification and pre-filtration, as well as sterile filtration. Top quality with high quality constructed housings. Teflon (PTFE) membrane material. Teflon (PTFE): hydrophobic membrane, chemically inert. Suitable for aggressive solvents, but also for air and gas filtration.

Housing material: polypropylene (PP)

Dead volume: 5 µl (4 mm), 10 µl (15 mm), 100 µl (25 mm)

Filter area: 0.07 cm<sup>2</sup> (4 mm).

## Additional Information

---

No.	8-7075
Manufacturer (Brand)	Ahlstrom
EAN	4058072051648
Transport temperature	Room temperature
Color	transparent (clear)
sterile	No
suitable for	Syringes
DM outside	4 mm
Filter properties	with special chemical resistance
Fluid behavior	hydrophobic
Material housing	Polypropylene (PP)
Material membrane	Polytetrafluoroethylene (PTFE. Teflon®)
Pore size	0.22 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer cone female
Area diaphragm	0.07 cm <sup>2</sup>
Dead volume	0 ml
Type filter	Syringe pre-filter
for medium	Liquids Organic solutions Aqueous solutions

