



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Ahlstrom ReliaPrep syringe filter, Teflon (PTFE),
D: 15 mm, non-sterile, pore size 0.2 µm, Luer
Slip, 50 pcs./pc**

**€181.80
plus VAT &
Shipping**

Product Images



Description

For sample preparation of aqueous and organic solutions, particle removal, high purification and pre-filtration, as well as sterile filtration. Top quality with high quality constructed housings. Teflon (PTFE) membrane material. Teflon (PTFE): hydrophobic membrane, chemically inert. Suitable for aggressive solvents, but also for air and gas filtration.

Housing material: polypropylene (PP)

Dead volume: 5 µl (4 mm), 10 µl (15 mm), 100 µl (25 mm)

Filter area: 1.7 cm² (15 mm).

Additional Information

No.	8-7078
Manufacturer (Brand)	Ahlstrom
EAN	4058072051679
Transport temperature	Room temperature
Color	transparent (clear)
sterile	No
suitable for	Syringes
DM outside	15 mm
Filter properties	with special chemical resistance
Fluid behavior	hydrophobic
Material housing	Polypropylene (PP)
Material membrane	Polytetrafluoroethylene (PTFE. Teflon®)
Pore size	0.22 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer cone female
Area diaphragm	1.7 cm ²
Dead volume	0.01 ml
Type filter	Syringe pre-filter
for medium	Liquids Organic solutions Aqueous solutions

