



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Ahlstrom ReliaPrep Syringe Filter, Regenerated Cellulose (RC), D: 25 mm, non-sterile, pore size 0.2 μ m, Luer Slip, 500 pcs./pc

€904.70
plus VAT &
Shipping

Product Images



Description

For sample preparation of aqueous and organic solutions, particle removal, high purification and pre-filtration, and sterile filtration. Top quality with high quality constructed housings. Regenerated cellulose (RC) membrane material. Regenerated Cellulose (RC): hydrophilic membrane with excellent chemical and solvent resistance. The filter for sample preparation in HPLC. Housing material: polypropylene (PP) Dead volume: 10 µl (15 mm), 100 µl (25 mm) Filter area: 1.7 cm² (15 mm) 4.8 cm² (25 mm).

Additional Information

No.	8-7070
Manufacturer (Brand)	Ahlstrom
EAN	4058072051600
Transport temperature	Room temperature
Color	transparent (clear)
sterile	No
suitable for	Syringes
DM outside	25 mm
Filter properties	with special chemical resistance suitable for HPLC
Area diaphragm	4.8 cm ²
Fluid behavior	hydrophilic
Material housing	Polypropylene (PP)
Material membrane	Cellulose. regenerated (RC)
Pore size	0.22 µm
Dead volume	0.1 ml
Type Connection Output	Luer cone male
Type Connection Input	Luer cone female
Type filter	Syringe pre-filter
for medium	Liquids Organic solutions Aqueous solutions

