



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Nalgene syringe filter w. CA membrane, 0.2 μ m,
25 mm \varnothing , sterile, 50 pcs/pack**

**€171.00
plus VAT &
Shipping**

Product Images



Description

Sterile disposable syringe filters

Diameter: 25 mm

Individually sterile packed, pyrogen-free.

Inlet side: Luer lock female Outlet side: Luer cone male.

For all standard syringes and cannulas.

PU: 50 pcs./pack

Additional Information

No.	CD-1877
Manufacturer (Brand)	Nalgene
old neolab article no.	3-1902
EAN	4058072219598
Transport temperature	Room temperature
Color	White
sterile	Yes
suitable for	Syringes
DM outside	25 mm
TBST MAX	121 °C
Filter properties	individually packed
MAX operating pressure	5 bar abs.
Fluid behavior	hydrophilic
Material housing	Polymethylmetacrylate (PMMA. PXM®) (acrylic glass)
Material membrane	Cellulose acetate (CA)
Pore size	0.20 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	2.8 cm ²
Dead volume	0.1 ml
Type filter	Syringe pre-filter
for medium	Organic solutions Aqueous solutions

