



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Syringe filter w. PTFE-Membran, 0.2 µm, 25 mm Ø, 100 pcs./pack

Special Price  
**€205.00** was  
**€307.50**  
plus VAT &  
Shipping

### Product Images



## Description

---

Syringe filter, housing made of glass fiber reinforced PP with PTFE membrane.

Can be used with aggressive chemicals including acids and non-aqueous solvents, e.g. GC and HPLC.

Also suitable for sterile aeration.

Non-sterile, sterilizable in autoclave or with gas.

## Additional Information

---

No.	3-1901
Manufacturer (Brand)	neolab
EAN	4058072018382
Transport temperature	Room temperature
Color	White
sterile	No
suitable for	Syringes
DM outside	25 mm
TBST MAX	121 °C
Filter properties	with special chemical resistance suitable for HPLC
Fluid behavior	hydrophobic
Material housing	Polypropylene (PP)
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.20 µm
Type Connection Output	Luer cone male
Area diaphragm	4.9 cm <sup>2</sup>
Type filter	Syringe pre-filter
for medium	Liquids

