



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**qpore® syringe prefilter, CA, glass fiber prefilter,
non-sterile, 0.45, Ø 30 mm, 100 pcs/pack**

**€129.00
plus VAT &
Shipping**

Product Images



Description

qpore® offers a comprehensive range of high quality **syringe filters** for various filtration applications in your laboratory. All filters are manufactured under the highest quality standards from the best raw materials.

This non-sterile qpore® syringe filter has a hydrophilic membrane made of CA and is therefore ideal for the filtration of aqueous solutions. CA exhibits low protein binding. Therefore, this syringe filter can be used wherever application or clear filtration under non-sterile conditions is allowed, such as mostly in sample preparation prior to HPLC or GC.

Prefilter made of glass fiber included.

The stable filter housing made of polypropylene is pressure resistant up to a maximum of 6.0 bar which allows fast filtration.

Features:

Low dead volume

Stable at pH 1-14

Luer connections: Luer lock female, Luer cone male

No risk of mix-up due to labeling (membrane type, pore size)

The syringe filters are packed non-sterile 100 pieces in a bag.

Additional Information

No.	6-0063
Manufacturer (Brand)	qpore
EAN	4058072171377
Transport temperature	Room temperature
Color	blue
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	30 mm
TBST MAX	50 °C
Filter properties	with glass fiber prefilter
MAX operating pressure	6 bar abs.
Area diaphragm	4.6 cm ²
Fluid behavior	hydrophilic
Material housing	Polypropylene (PP)
Material membrane	Polyethersulfone (PES)
Pore size	0.45 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Liquids

