




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



qpore® Syringe Filter, CME, sterile, 0.45 µm, Ø 33 mm

€179.00
plus VAT &
Shipping

Product Images



Description

This **sterile qpore® syringe filter** has a **hydrophilic membrane made of CME** with very good thermal and chemical resistance. With an excellent flow rate, it is optimally suited for sterile filtration, clear filtration and cell removal under sterile conditions. The effective filtration area of this syringe filter is 4.90 cm². The stable filter housing made of polypropylene is pressure-resistant up to a maximum of 5.0 bar, which enables rapid filtration.

Due to the Luer-Lock male connection variant, it is possible to screw several of these filters in series one after the other.

Features:

Membrane diameter 33 mm

Low dead volume

Stable at pH 1-14

DNA-, DNase-, RNase-, pyrogen-free

Luer connections: Luer lock female, Luer lock male

No risk of mix-up due to labeling (membrane type, pore size)

The syringe filters are individually sterile packed to 100 pieces.

Additional Information

No.	6-0077
Manufacturer (Brand)	qpore
EAN	4058072172602
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	Yes
suitable for	Syringes
DM outside	33 mm
TBST MAX	90 °C
Filter properties	suitable for HPLC single packed
MAX operating pressure	6 bar abs.
Area diaphragm	4.5 cm ²
Fluid behavior	hydrophilic
Material membrane	Mixed cellulose ester (MCE)
Pore size	0.45 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Liquids

