



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



qpore® Syringe Filter, CME, sterile, 0.45 µm, Ø 13 mm

€179.00
plus VAT &
Shipping

Product Images



Description

This **sterile qpore® syringe filter** has a hydrophilic **membrane made of CME** with very good thermal and chemical resistance. With an excellent flow rate, it is optimally suited for sterile filtration, clear filtration and cell removal under sterile conditions. The stable filter housing made of polypropylene is pressure resistant up to a maximum of 5.0 bar, which enables fast filtration.

The Luer-Lock male connection variant makes it possible to screw several of these filters in series one after the other.

Features:

Low dead volume

Stable at pH 1-14

DNA-, DNase-, RNase-, pyrogen-free

Luer connections: Luer lock female, Luer lock male

No risk of mix-up due to labeling (membrane type, pore size)

The syringe filters are individually sterile packed to 100 pieces.

Additional Information

No.	6-0128
Manufacturer (Brand)	qpore
EAN	4058072173227
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	Yes
suitable for	Syringes
DM outside	13 mm
TBST MAX	90 °C
Filter properties	suitable for HPLC single packed
MAX operating pressure	6 bar abs.
Fluid behavior	hydrophilic
Material membrane	Mixed cellulose ester (MCE)
Pore size	0.45 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	1.09 cm ²
Type filter	Syringe pre-filter
for medium	Liquids

