




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



qpore® Syringe Filter, CME, non-sterile, 3.0, Ø 13 mm

€99.00
plus VAT &
Shipping

Product Images



Description

This non-sterile **qpore® syringe filter** has a hydrophilic membrane made of CME with very good thermal and chemical resistance. With an excellent flow rate, it is optimally suited for clear filtration.

- Filter housing made of polypropylene, maximum pressure resistance 5.0 bar
- Fast filtration

Additional Information

No.	6-0137
Manufacturer (Brand)	qpore
EAN	4058072175450
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	13 mm
TBST MAX	90 °C
MAX operating pressure	6 bar abs.
Fluid behavior	hydrophilic
Material membrane	Mixed cellulose ester (MCE)
Pore size	3 µm
Type Connection Output	Luer lock male
Type Connection Input	Luer lock female
Area diaphragm	1.09 cm ²
Type filter	Syringe pre-filter
for medium	Liquids

