



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## qpore® Syringe Filter, CME, non-sterile, 8.0, Ø 13 mm

**€99.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

This non-sterile **qpore® syringe filter** has a hydrophilic membrane made of CME with very good thermal and chemical resistance. With an excellent flow rate, it is optimally suited for clear filtration.

- Filter housing made of polypropylene, maximum pressure resistance 5.0 bar
- Fast filtration

## Additional Information

---

No.	6-0139
Manufacturer (Brand)	qpore
EAN	4058072175474
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	13 mm
TBST MAX	90 °C
MAX operating pressure	6 bar abs.
Area diaphragm	1.09 cm <sup>2</sup>
Fluid behavior	hydrophilic
Material membrane	Mixed cellulose ester (MCE)
Pore size	8 µm
Type Connection Output	Luer lock male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Liquids

